Release date: 27.9.2019
No.: 0200151

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Trdota po Vickersu: > 26 HV
Upogibna trdnost: > 114 MPa
E-Modul: > 2771 MPa
Prosta monomera: < 1%

DIMENZIJE

Premer: 95, 98, 98.5, 100 mm
Debelina:

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DIMENSIONS

Diameter: 95, 98, 98.5, 100 mm
Thickness:

- Basic, Splint: 12, 13, 14, 15, 16, 18, 20, 25, 30 mm
- Multilayer: 15, 16, 18, 20 mm
- 3-layer: 15, 16, 18, 20, 25 mm

POLIDENT PMMA CAD/CAM DISC TYPE BASIC (MONOCHROME), MULTILAYER AND 3-LAYER

POLIDENT PMMA CAD/CAM DISC type BASIC (monochrome), MULTILAYER and 3-LAYER are suitable for manufacturing of long-term temporary crowns and bridges and for verifying the bridge on the plaster model/mouth before the final ZrO2 milling process.

POLIDENT PMMA CAD/CAM DISC type MULTILAYER and 3-LAYER are composed by several different layers. Layers colours are adapted to the natural colour gradient of the tooth. The final result is a very aesthetic and economic restoration.

SOME RULES TO CONSIDER THROUGH THE PMMA C & B DESIGN:

- The minimum wall thickness has to be 1,0mm in occlusal region and 0,6mm in cervical region.
- The connector profile has to be more than 10mm² in anterior region and more than 15 mm² in posterior region.
- A max of two pontics is allowed between two abutment teeth in posterior region.
- Pay attention to the positioning of the restoration within the MULTILAYER and 3-LAYER disc. The enamel layer is marked on the top of disc.

FINISHING C & B AND CEMENTING IN THE MOUTH:

- Use burs for acrylic. Polishing is possible with polishing paste and soft hair brushes. Pay attention not to overheat the material!
- For cementing the bridge use temporary cement.

SHADES

- Basic: A1 / A2 / A3 / B1 / B2 / B3 / BL1 / BL2 / BL3 / E1 / E2
- Multilayer: A1 / A2 / A3 / B1
- 3-layer: A1 / A2 / A3 / B1

POLIDENT PMMA CAD/CAM DISC TYPE SPLINT

POLIDENT PMMA CAD/CAM DISC type SPLINT is mainly suitable for fabrication of splints and other oral cavity devices. Discs can also be used for fabrication of frameworks for testing within the oral cavity, for working models and as a base for fabrication of moulds in precision casting with metal or press ceramics. Braces are first designed with software for computer-aided design for braces. Thickness and cross-sections that apply to the final material need to be taken into consideration during design of frameworks for precision casting.

Tools for acrylate are used for finishing and final adjustments. For polishing, a pumice stone, pastes and soft polishing brushes are used.

It is important not to overheat the material!

SHADES

- Splint: Clear / Berry / Blue / Green / Orange / DE

STORAGE

Store at room temperature.

WARNING

In case of known allergies or hypersensitivity to acrylic materials do not use the product.

In case that any serious incident has occurred in relation to the device it should be reported to the manufacturer (Polident d.o.o.; vigilanca@polident.si) and the competent authority of the Member State, where the user is established.

CE 2460

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Vickershärte: > 26 HV
Biegefestigkeit: > 114 MPa
E-Modul: > 2771 MPa
Restmonomer: < 1%

GRÖSSE

Durchmesser: 95, 98, 98.5, 100 mm
Dicke:

- Basic, Splint: 12, 13, 14, 15, 16, 18, 20, 25, 30 mm
- Multilayer: 15, 16, 18, 20 mm
- 3-layer: 15, 16, 18, 20, 25 mm

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